Multi-lingual information for parents on how to use eye patches

to successfully treat visual impairments

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The most common causes of visual impairment are:

- Short- and long-sightedness
- Squinting
- Corneal curvature
- Different types of defective vision affecting both eyes
- Ocular tremor

What are possible causes of visual impairment?

Once a child has passed through childhood the sense of sight is fully developed.

If one or both of your child's eyes do not reach full visual acuity in spite of them wearing corrective glasses, this is referred to as 'amblyopia'. Full development of the sense of sight can be impaired by squinting. Children suppress the squinting eye's visual impression, consequently visual acuity cannot develop in the weaker eye.

How can visual impairments be treated?

Visual impairments can be treated during childhood using eye patch therapy and glasses. To train the weaker eye, the eye that sees properly is covered with an eye patch. The ophthalmologist determines length of the therapy, as well as how often and how long the eye needs to be covered, on an individual basis for each child. If the child wears glasses, they must be worn at the same time as the eye patch! Beneficial activities while the child is wearing the eye patch could be: painting, making crafts, writing, reading, solving puzzles, playing board games.

What are the possible consequences of untreated visual impairment?

If eye patch therapy is not carried out consistently in their childhood, the visual impairment will persist throughout their entire life. The visual impairment **cannot** be treated later in life by laser therapy, surgery or any other type of treatment. An untreated visual impairment may restrict the range of professions a child would like to pursue in adulthood.

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

Instructions and tips for applying eye patches

How do you apply the eye patches?

Please ensure that your **child's skin is clean, dry and free of any creams or salves**. Pull the white paper film off the eye patch.

1. Place the patch on the **closed eye** such that the **white inner lining is over the eye socket**.
2. Position the **nose slot centrally over the nose**.
Now press down the edges of the patch so they stick.
3. If your child has **glasses**, these should always be **worn over the eye patch**.
4. After the recommended time, gently **remove the patch working from the nose** to the edge of your child's face.


Tips:

Children have different skin types and the eye patches may stick better to some types than others.

Piratoplast® eye patches are coated with a **hypo-allergenic adhesive** and have been **positively tested to be gentle on the skin**. Despite this and depending on the time of year and what your child likes to get up to, patches may stick to varying degrees, e.g. if your child sweats a lot during the summer or when romping around.

If you believe that the **eye patch is sticking too well**:

- briefly first apply it to **back of your hands** before sticking it over your child's eye

I put an eye patch over my left eye so that I can train the right one to see better.

Discover more

You will find everything you need to know about eye patches at:
www.piratoplast.de/ratgeber

- Useful information for parents about treating eyes with eye patches
- Accompanying materials for children (e.g. instructions for handicrafts, collection poster, activity books and much more)

Have a look at our Facebook page:

**[www.facebook.com/
Piratoplast.Augenpflaster](https://www.facebook.com/Piratoplast.Augenpflaster)**

Discover our digital collective poster and visualise the successes of your little pirate
www.piratoplast.de/login